IMPLANT

Neo 

TÉCNICA QUIRÚRGICA

Clavo placa pediátrico

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
Correction Angle Calculation: Calculate the correction angle (varus) by subtracting the current neck-shaft angle (measured via fluoroscopy in 30-40° internal rotation) from the desired angle. Add this result to the selected plate angle (90° or 100°) to obtain the final guide wire insertion angle.

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Guide Wire Insertion: Attach the angle guide to the lateral femur with a clamp. Advance the guide wire 10–15 mm and verify optimum positioning via AP and lateral fluoroscopy to avoid cortical penetration. Once confirmed, advance the wire to within 5–10 mm of the calcar or proximal physis.

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Insert the selected chisel over the guide wire using the chisel guide. Verify proper sizing and centered alignment through the femoral neck isthmus via AP and lateral fluoroscopy. Drive the chisel with a mallet to the planned depth (5–10 mm short of the physis or calcar), ensuring the wire does not bend or enter the joint. Remove the chisel guide and back out the chisel 5–10 mm to facilitate later removal.

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


Use the osteotomy gauge to mark the initial cut level above the lesser trochanter, verifying the line via fluoroscopy. Advance the oscillating saw slowly to make the first transverse cut, controlling the proximal fragment with the chisel. For precise derotation, make a vertical reference score on the anterolateral cortex beforehand. If required, make a second transverse cut for femoral shortening, and finally remove the seating chisel.

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Select the largest and most stable implant that fits the femur to allow early mobilization without a cast. Balance distal fragment medialization based on varus correction (higher in CP than DDH/LCP). Attach the plate to the threaded inserter and advance it over the guide wire. Ensure the wire is not bent to prevent the implant from locking or driving it into the joint.

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Attach the implant to the threaded inserter and advance it over the guide wire. It is mandatory to check that the guide wire has not bent during previous steps to prevent instrument locking or accidental displacement.

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Reduce the osteotomy by manipulating the distal fragment (flexing it to align with the proximal fragment). Secure the plate against the lateral femur using a bone clamp in the second-to-last shaft hole. Verify via fluoroscopy that the plate is parallel and centered on the femoral shaft to prevent flexion or extension deformities and ensure bicortical fixation.

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
Place the gold end of the drill guide into the top compression hole with its arrow facing the osteotomy. Drill through both cortices and use the depth gauge to hook the far cortex for screw measurement. Always insert compression (non-locking) screws first, followed by locking screws, leaving the proximal locking screw for last.

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
Attach the hexalobe driver to the in-line ratchet in forward position. Insert the selected cortical screw through both cortices. Note that due to dynamic compression, the screw may protrude 1–2 mm more than expected, which is acceptable. Repeat steps 8 and 9 for the second compression hole.

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Remove the bone clamp. Insert the threaded drill guide into the locking plate hole, avoiding cross-threading. Insert the calibrated drill bit and drill through both cortices. Read the screw measurement from the calibrated drill bit off the back of the threaded drill guide.

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Insert the measured locking screw using the hexalobe driver and ratchet. Repeat for the second shaft screw, and place the proximal locking screw last. Drill the proximal screw hole under mandatory fluoroscopic control to avoid femoral neck perforation; this screw runs parallel to the blade and is typically 10–15 mm shorter.

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Following final fixation, assess the stability of the implant, the reduction, and the hip joint both clinically and via intraoperative fluoroscopy (an arthrogram may be used). Finally, irrigate the incision and perform a layered closure.